DEUTSCH

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ENGLISH

Important safety information – retain for later use.

In order to make optimum use of all the merits of the clinical thermometer, you should carefully read the operating instructions before using it, retain the instructions for later use and also keep them accessible for other users.

- The thermometer is intended to measure the temperature of the human body.
- Do not allow children to play with the device. Children must be supervised at all times whilst using the device.
- The comforter thermometer is intended only for oral measurement of body temperature in children up to 5 years of age. It is not suitable for constant sucking.
- Before every use, check the comforter thermometer and especially the silicon nipple for signs of damage or wear. Damaged or worn comforter thermometers should not be used.
- Check that the teats are firmly attached to the pacifier. There is a choking hazard.
- Contact your doctor immediately if any components are swallowed.
- Under specific conditions, the displayed temperature may be considerably different from the actual body temperature of the child (-0.5°C to 0.5°C). Therefore, please ensure that the child does not eat or drink anything shortly before the measurement, and that the child is not crying or breathing through the mouth when the temperature is being measured.
- In case of doubt and signs of fever (temperature higher than 38.5°C, flushed appearance, heat, perspiration, restlessness) measure the body temperature afterwards by using a conventional clinical thermometer in the anus (rectal) and discuss the measured temperatures with the physician.
- Always wait for the acoustic signal which indicates the minimum time required for measurement. Discuss the temperature with your family physician.

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- The thermometer contains sensitive electronic parts. Therefore, protect it from blows, bending, high temperature or direct exposure to sunlight.
- Keep the thermometer away from solvents and other chemicals.
- When the thermometer is switched on, it carries out a self-test. Therefore the accuracy of measurement need not be tested at regular intervals.
- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. For more details, please contact customer service at the address indicated.

Use

The pacifier is delivered in one size. Before using it for the first time, check that the pacifier fits into the mouth of your child. If the pacifier is too large, there is a choking hazard. Press the ON switch to turn the unit on. A short acoustic signal will announce that it is switched on. The thermometer will carry out a self-test for about 2 seconds. During the self-test all segments of the display will be visible. Then "Lo °C" or the value of the last measurement will appear in the display. Make sure the baby keeps the comforter in the mouth during the measurement and breathes calmly through the nose. The current temperature will be shown directly, e.g. 32.1°C. The measurement takes approximately 5 minutes, which is indicated by an acoustic signal. By prolonging the length of measurement until after the acoustic signal, more accurate measurements can be obtained. If temperatures higher than 42.9°C are measured, the symbol for too high temperature, "Hi" will appear. In order to increase the life of the battery, please switch off the thermometer by briefly pressing the ON key after the measurement.

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